# **Operation Manual**

# Isotemp RT AVCD HPS Advanced Hotplate Stirrer



9240-11-012 V01 12/01/17



# This manual cover the model is shown below

NA Model	EU Model	Voltage	Description
11676261	15326607	230V-EU, UK, ANZ/CN	Isotemp RT AVCD HPS230V
11676262	N/A	120V-US	Isotemp RT AVCD HPS120V



Important Before using this product, read this entire operation manual carefully. Users should follow all of the operational guidelines contained in this manual and take all necessary safety precautions while using this product. Failure to follow these guidelines could result in potentially irreparable bodily harm and/or property damage.

Caution All internal adjustments and maintenance must be performed by qualified service personnel.

Material in this manual is for information purposes only. Fisher Scientific is committed to a continuing program of product development and improvement, and reserves the right to change information, such as specifications, appearance, and dimensions, described in this document without notice. Fisher makes no representations or warranties with respect to this manual. In no event shall Fisher be held liable for any damages, direct or incidental, arising out of or related to the use of this manual.

No part of this manual may be reproduced or transmitted in any form or by any means, including photocopying, recording, or using information storage and retrieval systems, for any purpose other than the purchaser's own use, without the express written permission of the manufacturer.

Any other product names and services identified in this manual are trademarks or registered trademarks of their respective owners. No such use, or the use of any trade name, is intended to convey endorsement or other affiliation with Fisher Scientific.

©2017 Fisher Scientific. All rights reserved.

Preface

This manual contains important safety and operation information. You must carefully read, understand, and follow all the instructions in this manual prior to operating this instrument. Keep this manual in a safe place nearby for reference and make it easily available to all users.

- 1) This manual highlights DANGER/WARNING/CAUTION/NOTICE alerts to prevent injury or property damage and also to achieve optimum performance of your instrument.
- 2) These alerts are classified into four types in this manual depending on the importance and the risk levels as described below:

Symbols	Meaning
	Indicates a hazardous situation which, if not avoided, will result in death or serious injury.
<b>A</b> WARNING	Ignoring this warning could cause serious injury or even death.
	Ignoring this caution could cause injury or property damage.
NOTICE	Ignoring this notice could cause operational problems.

- 3) The claim which is out of the quality guarantee published by the Manufacturer is out of Manufacturer's responsibility.
- 4) The damage which is from unexpected fault or damage of user by Acts of God is out of Manufacturer's responsibility.

# **Table of Contents**

Section 1	Warning and Cautions	1-1
Section 2	Functional Description	2-1
	Features	2-1
	Safety	2-2
	Convenience	2-3
	Construction	2-5
Section 3	Unpacking and Installation	3-1
	Location Conditions	3-1
	Pre-startup Check	3-2
	Connecting to Power Supply	3-2
Section 4	Control Panel	4-1
	Check Power Status	4-3
	Connecting the Temperature Probe	4-3
	Stirrer Knob	4-5
	Run/Stop Button	4-5
	Temperature Setting	4-6
	Temperature Setting During Operation	4-7
	Resetting Set Temp During Operation	4-8
	Maximum Heating Rate Setting	4-9
	Timer Setting	4-10
	Timer Setting During Heating	4-11
	Check Remaining Run Time	4-12
	Temperature Probe Control Mode	4-12
	Configuration Mode	4-14
	Temperature Unit Conversion( $^{\circ}C \leftrightarrow ^{\circ}F$ )	4-14
	Selection of the Timer Mode	4-16
	Temperature Offsetting(OFST/OFSP*)	4-17
	Maximum Heating Rate Setting(LIMt)	4-19
	Change the High Temp Limits(t-H)	4-20

	Change the Low Temp Limits(t-L)	4-22
	Check Configuration Settings(LISt)	4-23
	Auto-tuning(PID Parameter Calibration)	4-24
	Coefficient(COEF)	4-25
	Default(DEFA)	4-27
	Motor Minimum(Mmin)	4-28
	Escape Mode(ESC)	4-29
Section 5	Safety Device	5-1
Section 6	Maintenance	6-1
	Cleaning Product	6-1
	Relocation	6-2
	Keeping Product	6-2
Section 7	Trouble Shooting	7-1
	Error Messages	7-3
Section 8	Accessories	8-1
	Assembly of Support Rod	8-1
	Transparent Shield & Heating Bath Assembly	8-2
Section 9	Technical Specifications	9-1
	Disposing of the Unit	9-2
Section 10	Warranty Information	10-1

# Section 1 Warnings and Cautions

#### **A WARNING**

# Ignoring the following warnings could cause serious injuries or even fatal accidents.

Check and connect properly -the voltage, phase and capacity of power supply on the ID plate before installation.

Power supply must be properly grounded. Abnormal grounded connection causes serious damage. Grounded connection must not be on the water pipe and gas pipe.

Use correct and provided power code

Do not install the product in the place that the gas could leak out. Do not use in the place that has industrial oil smoke and/or metallic dust. It causes fire or electric shock.

Do not use the machine near to places where explosion can occur due to organic evaporating gases.

Explosive materials: Acid, Esther, Nitro compound.

Inflammable materials: salt peroxides, inorganic peroxide, salt acids.

Check equipment for permissible environmental condition. It can be the cause fire or trouble by electricity, electronic, and damage of motor.

Permissible environmental condition - Temperature 2°C to 60°C, Maximum relative humidity 80%.

Wear your personal protective equipment in accordance with the hazard category of the medium to be processed; splashing liquids, projectile parts, body parts, hair, clothing and jewelry getting caught.

Unplug, if there is strange sound, smell and smoke from the product. Stop operating and request service.

Keep out of direct sunlight. It may influence product life and proper operation.

Do not use the machine in places where moisture is high and flooding can occur.

Do not assemble, repair, modify on your own. The product may not work well and electric shock in the efficiency of the product. Also you will void the warranty.

#### 

# Ignoring the following cautions could cause injuries or property damages.

Indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.

Do not touch the top plate or any object near it even when the heater is turned off. You may get burned due to residual heat.

After using this equipment, make sure to turn off the main power switch and also to disconnect the power cord from the power outlet for the safety of other users.

Do not place heavy objects, including this equipment, on top of the power cord and do not strip, scratch, bend, twist, pull, or heat the power cord. A damaged power cord is a fire and electrical shock hazard.

Make sure to set up this equipment on a flat, stable, clean, non-slip, dry, and fireproof surface inside a lab with proper safety measures.

Do not place any device which can be affected by the motor vibrations near this equipment.

Do not touch the power outlet, power socket, or power cord with wet hands. And make sure to connect the power cord directly and firmly to the power outlet and power socket.

Do not put or insert any objects (especially if conductive or flammable) inside this equipment.

Do not expose this equipment to any heat sources including direct sunlight.

Beware that mechanical shock or vibration can damage this equipment and pay extra attention while moving it. Damages caused by mechanical shock or vibration may result in injury or fire.

Do not impact the top plate or heat sink. You can damage the equipment or get injured.

Do not install this equipment near any device that generates high frequency noise such as high frequency welding machines, high frequency sewing machines, or SCR power controllers.

Before cleaning, make sure to unplug the power cord to avoid electric shock or fire.

Do not use chlorine bleach, ammonia-based cleaners, abrasives, ammonia, or metal scouring pads. Wipe with a soft damped cloth or a sponge soaked in water or diluted neutral detergent.

The Fisher Scientific Isotemp RT AVCD HPS model offers temperature and control by external temperature probe. This unit also boasts quick heat-up time thanks to ample heating capacity as well as the superb heat transfer rate enabled by the tightly integrated structure of the heater and the ceramic-coated aluminum alloy top plate. Also provides a variety of temperature modes relating to temperature control, heat-up time.

#### Features

- By adopting a special magnet with exceptionally strong magnetic coupling power as well as smooth-start stirring mechanism. You will experience virtually no decoupling of magnetic stir bars even with viscous media or at high speeds (30~2000rpm, guaranteed). The BLDC motor and the special magnet also provide big stirring power.
  - · Real-time stirring with quick operation or stop
  - When the external temperature probe is coupled, it automatically changes to external sensor mode. User can check media temperature also control.
  - Microprocessor PID Feedback Control

Fast and precise temperature control is provided by the microprocessor PID controller.

Selection of the Temperature Control Modes

Three user-selectable temperature control modes are provided for your convenience; Optimum, Fast, and Slow.

• Auto-Tuning

Automatic tuning of the PID parameters provides more accurate temperature control.

Changing Temperature Limit

The top plate temperature control range can be changed to protect media by changing over temperature limit value or low temperature limit value.

Features	<ul> <li>Offset To use your own thermometer for temperature control in specific applications, there can be some differences between the temperature of your thermometer and the displayed temperature of this unit. If needed, you can offset such temperature differences ranging from - 10°C to ±50°C at 0.1°C intervals. </li> <li>Selection of the Timer Mode Immediate Activation (T1): The timer starts immediately after setting the timer, Delayed Activation (T2): The timer is activated only when the set temperature is reached. </li> <li>Quick Heat Up Time Ample heating capacity and superb heat transfer rate of the tightly integrated structure of heater and ceramic-coated aluminum alloy top plate allow quick heat up time. The ceramic-coated top plate is highly resistant to heat and corrosion. In addition, its white color is optimal for monitoring color changes of the media during operation.</li></ul>
Safety	<ul> <li>Hot Top Warning Indicator The top plate temperature can remain very hot for some time even after the heater is turned off. To prevent injury or fire under such circumstances, this instrument has a hot top warning indicator on the control panel. This indicator will illuminate if the top plate temperature is over 50°C. Even so, do not rely on this indicator alone for your safety </li> <li>Overheating prevention device Basically built in electronic overheating prevent device and mechanical body overheating prevent device. </li> <li>Heating plate overheating prevention device A heating plate overheating prevention device is particularly built-in. </li> <li>The devices WARNING LED will be lighted. In case of the devices temperature goes over safe range.</li> </ul>

Safety


Error indication regarding external temperature sensor

In case of Error message, the Probe does not contain enough in media or on air.

- A transparent shield is also provided to allow you to monitor your operations more safely by shielding you against decoupled stir bars or liquid splashes. (optional)
- Heating Bath (optional)

A specially designed non-slip heating bath comes with the instrument. To prevent unintentional slips of the bath during operation, this bath has concavo-convex bottom shown below so that it can sit stably on top of the top plate.





# • Conveniently set indicating temperature and stirring speed by two different knobs:

- In case of external temperature probe coupled, automatically change to external sensor mode.
- The external temperature probe (sensor) is made by Pt material. So, the probe (sensor) can be fix easily on lab ware.
- Showed on digital display- top plate temperature, Probe measuring temperature, stirring speed, timer setting value, and remaining time value.
- In case of the stirring media viscosity change, the feedback function maintains regular speed.

## Convenience

## Convenience

- Smooth-start stirring system reduces magnetic bar decoupling.
- Two support rods are provided to hold various kinds of devices such as temperature sensors, thermometers, laboratory glassware, etc. (optional).
- User designation temperature range

Your own temperature limit can be easily set and protect simple from controlling time or user mistake.

- User designation stirring speed range
  - User can set up the lowest stirring speed. So the device can be stirring as quickly as possible time.

## Construction




Figure 2-2. Front Components

(1) Temperature probe

This is PT100 sensor to check media temperature.

(optional 2 different temperature probe series. Refer to Accessories list)

- (2) Heat sink
- (3) Heater knob
- (4) Control panel
- (5) Stirrer knob
- (6) Clamp holder (optional)
- (7) 3 Prong clamp (optional)
- (8) Support rod (optional)
- (9) Heating bath (optional)

## Construction




Figure 2-3. Back Components

- (10) Power switch
- (11) Transparent Shield (optional)
- (12) Top plate
- (13) Threaded hole
- (14) Power socket
- (15) Temperature probe cap

# Section 3 Unpacking and Installation

Upon receiving the instrument, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred during shipment must be detected before unpacking. If such damage is found, notify the carrier immediately.

After unpacking, check to ensure that all the following parts and accessories are included in the package. If not, contact your dealer or Fisher Scientific immediately.

ltem	Figure	Quantity	Description
Main body		1	
Power cord		1	
Operation Manual	A.	1	
Magnetic bar	Ø	1	-
Temperature	() <sup>r</sup> O	1	B class
probe	are a	I	Max 250°C

## **Location Conditions**

Place the unit on the wrench, or table when in use, and observe minimum distances in 30cm between the other devices.

#### 

The unit should be located away from naked flame sources, direct sunlight. It can come with the malfunction or lower the function.

Section 3 Unpacking and Installation

### **Pre-startup Check**

- Check the unit is balanced well.
- Never install or use this instrument with or near to hazardous or flammable substances.
- Do not install this instrument near any device that generates high frequency noise.
- Do not install with short circuit, water leak, and risk of flooding places.
- Do not install with industrial harmful gas, metal dirt environments.

# Connecting to Power Supply

When connecting power, use only the power cord that came with your instrument. The power connection procedures are as follows:

- Before connecting the power cord, make sure that the main power switch is turned off.
- (2) Plug the power cord into the power socket at the back of your instrument as shown in the diagram below.
- (3) Plug the other end into a properly grounded and dedicated power outlet nearby.




Figure 3-1. Connect Power Cord


#### **A WARNING**

**Electrical Shock Hazard.** 



- Check voltage, phase and capacity of power supply and connect properly.
- Do not insert damaged line cord or multiple plugs into outlet at the same time.
- DO NOT use without safety PPE (working clothes, gloves, glasses)
- DO NOT handle or touch electrical cord or electrical parts with wet hands.
- Make sure to connect this instrument only to properly grounded power outlets to protect you and your instrument.

**Control Panel** 



Section 4 Control Panel

# **Control Panel**


		When this o	device is stop conduction, Configuration		
		Functions can be set. When user wants to change			
		mode, press the Config. Button with stirrer knob to			
		clockwise turning. Then the mode can be switching in			
		the following order.			
		Unit	Changing the temperature unit "Cels" -		
		Unit	Celsius, "Fahr" - Fahrenheit		
		Time	Changing timer mode		
		Offset	Changing Offset setting		
		11	Set maximum heater rate when used to		
		Limit	temperature control.		
			Highest temperature set top plate and		
0	CONFIG	t-H	external temperature sensor		
G	button	t-L	Lowest temperature set top plate and		
			external temperature sensor		
		List	Indicating the set value to turning stirrer		
			knob.		
		Calibration	Matching Auto tuning the device to		
			current situation or best condition.		
		Coefficient	Select to using control temperature		
			coefficient.		
			Change all the settings to the default		
	Default	values			
		Motor	Minimum stirring speed set (range :		
		Minimum	1~1000rpm)		
		Escape	Config. mode to be end		
Н	Speed	l display	STIRRER speed, timer		

Knob – Set a value and select a function, turn knob left or right.

Ι	Heater knob	Temperature setting
L Stirror knob	Stirror knob	STIRRER speed setting, value setting on
J	J Stirrer knob	Config.

### **Check Power Status**


(1) The following displays will appear in sequence together with a beep sound when the power cord is connected:



• Note that the following display will also appear momentarily if the main power switch is off:




(2) If you turn on the main power switch, the following display will appear in sequence:




# Connecting the Temperature Probe

When the external temperature probe connects with body, the temperature display part shows Probe measuring temperature. If the probe disconnect with body, the temperature display part shows top plate temperature.

- (1) Open the external temperature probe cap and turn right.
- (2) Insert external temperature probe with body.




(3) When you connect with temperature probe, the LED lamp is automatically turned on.




Section 4 Control Panel

# Connecting the Temperature Probe

Temperature display windows are showed on "Sensor change" message. After then, the display windows showed "Probe" measuring temperature.



(4) In case of separating temperature probe from main body, the temperature display windows are showed follow message;







- When you connect with temperature probe, check first hotplate temperature.
- When you connect with temperature probe and hotplate temperature is over 50°C, the "HOT TOP" LED will be lit.




- Pay attention to a high temperature.
- Set the Probe under 20mm from Media. If Probe is not set perfectly or on air, the Probe will be "Error" or the hotplate will be overheating.

#### Stirrer Knob

When the stirrer knob is turned right, the device will start operating. And the stirrer function is operating individually without heating part.

[Turning STIRRER operation, State indication LED]



During STIRRER function operating, the stirrer value will be showed on display window. When you push the stirrer knob, you can check current stirrer operating value.

In case user wants to stop stirrer operation, just turn knob all the way left.

#### NOTICE


 Stirrer speed range 1~2000rpm, you can control more correctly from 30 to 2000rpm.

# Run/Stop ButtonEven if the main power switch is turned on, you have to press the<br/>RUN/STOP Button to start or to finish the operation.

The RUN LED illuminates only when the unit is operating and the HEAT LED blinks only when the heater is activated.

If you press the RUN/STOP Button to start the operation and the current temperature is below the set temperature, both the HEAT and the RUN LEDs will light up as shown below.

[During heater operating, the state indication LED]



Note, however, that the HEATER LED will be turned on and off during the operation because the heater is automatically activated or deactivated to maintain the set temperature. Section 4 Control Panel

## **Run/Stop Button**





#### **A CAUTION**

- Hot Top Warning Indicator will illuminate if the top plate temperature is over 50°C.
- Be aware and careful of surface temperature.
- Even if the instrument is turned off, the surface of the top plate and the vessel on top of it will remain very hot for some time. Never leave your instrument accessible to others while it is hot and never touch it unless you are absolutely sure.

### NOTICE

- You can check set temperature during heater operating. Press the heater knob.
- In case of user requirements, press TEMP button to change set temperature during operating.
- User can check the elapsed time during operation by pressing TIMER button.
- If user wants to go TIMER setting mode during operation, just press timer button one more time.


#### **Temperature Setting**

The temperature setting procedures are proceeded by using TEMP button and the Heater Knob (left control knob).

Step 1: Press the TEMP Button and check the display showing the current temperature setting.



Step 2: If you want to change the temperature setting, select the desired temperature by turning and press the Heater Knob.



4-6

Section 4 Control Panel

#### **Temperature Setting**

Step 3: Push the RUN/STOP Button.

Heater LED will be on and heater will be activated.



If the external temperature probe is connected, the temperature which probe detects will be displayed. If not, top plate's temperature will be displayed.


### NOTICE

- Changing the temperature setting is allowed only within the low and high temperature limits. If changing the temperature setting cannot be done properly, check the low and high limits first. (Refer to Check the Configuration Settings)
- You can set the low and high temperature limits at any value between 0°C and 350°C.

# Temperature Setting During Operation


Press TEMP Button or the Heater Knob during operation to see the set temperature.

If you press TEMP Button, set temperature is shown on the right. After 15 seconds, the mode will be escaped automatically.



If you press the Heater Knob, set temperature is shown on the left. After


3 seconds, the mode will be escaped automatically.




Section 4 Control Panel

# Resetting Set Temp During Operation

Step 1: Press TEMP Button during operation.




Step 2: Select the desired temperature by turning and press the Heater Knob.



After 15 seconds, the mode will be escaped automatically.

Or push RUN/STOP Button as shown below to escape the mode.

Step 1: Push RUN/STOP Button during resetting the set temperature.



Step 2: Turn the Heater Knob clockwise and press the Heater Knob as shown below.



Section 4 Control Panel

# Maximum Heating Rate Setting


The higher the heating rate, the faster the heat up time but the wider the temperature overshoot and undershoot. Therefore, if you want to reduce the temperature fluctuation, you need to limit the heating rate to a certain degree by limiting maximum heating rate. [Same function as 4-19]

If, for example, the current heating rate is 100% but you want to set the heating rate limit at 50%, do as follows:


Step 1: Press the CONFIG Button. The default value is 100%.




Step 2: Change the limit to the desired value by turning the Stirrer Knob and press the Stirrer Knob.



Step 3: Save the desired limit by using the Stirrer Knob as shown below. [When you need to save the changes]



If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]



Step 4: Terminate the configuration mode as described in Configuration Mode.

### **Timer Setting**

This unit provides two types of timer mode: immediate timer activation and delayed timer activation. (For selecting the timer mode, Refer to 4-16 Selection of the Timer Mode)

[Indicating LED status during timer is activating]



You can stop timer by pushing RUN/STOP Button during timer's activation. For both timer modes, the timer setting procedures are the same:

Step 1: Press the Timer Button to begin the timer setting and also to check the display showing the current timer setting. Note that the hour frame is blinking first as shown below.




Step 2: If you want to change the timer setting, select and then save the desired values into the hour frame and the minute frame by turning the Stirrer Knob and pushing it for confirmation as shown below.



Step 3: Now, the unit will display the set temperature for verification as shown below. Turn the Stirrer Knob to select the desired value and then press it for confirmation. The following display will appear momentarily.


### **Timer Setting**

#### [If you change the set temperature]




Step 4: Then, the following display will appear asking whether you want to start the timer operation or not. If you want to save the new timer setting and to start the timer operation, you are required to change 'NO' to 'YES' by turning the Stirrer Knob and press it for confirmation as shown below.

[Operating confirmation]



Step 5: When the timer operation ends, you will be alerted by audible signals as well as the following display.



Terminate the timer operation either by pressing the RUN/STOP Button or by turning off the main power switch.

# Timer Setting During Heating

Press TIMER Button once to check run (remained) time of heating if heater is in operation.



Press TIMER Button twice to set timer during heating and follow above instruction (1) for timer setting.

# Check Remaining Run Time

Check the remaining run time during a timed operation by performing the following.

(1) Press TIMER Button once to check the remaining time during timed operation.




After 15 seconds, the mode will be escaped automatically.

Escape of timer mode is shown as below.

Step 1: Press RUN/STOP Button.


ESE	ם ה
-----	-----

Step 2: Turn the Stirrer Knob clockwise and press the Stirrer Knob.



Temperature Probe Control Mode Connect temperature probe to the unit. The default mode will be automatically changed to the Temperature Probe Control Mode and display shows the temperature value detected by the probe.

[For the temperature probe connection to the unit, Refer to Connecting Temperature Probe]



4-12

Section 4 Control Panel

# Temperature Probe Control Mode

For temperature and timer setting on the Temperature Probe Control Mode, see Temperature Setting and Timer Setting.

In case of low smoke point media heating, High temperature limit value should be lower than smoke point. [Refer to High Temperature Limit Setting]

### (1) Temperature Probe Errors

On temperature control mode if the external probe is in the air or is not enough immersed into the media, error message will be shown on display and heater will be automatically turned off.

For correct control media by the external temperature probe, make sure the sensor to completely immerse into the media (more than 20mm deep). If the sensor is in the air due to evaporation of media, precise temperature control may fails.

Error code	Cause		Effect
ERRO NO 6		If the temperature probe	Heater off
	Temperature Probe Errors	is in the air	rieater on
ERRO NO 7		If the temperature probe	
		is not completely immer	Heater off
		sed into the media	

#### [Refer to Error Message]

### 



- Hot Top Warning Indicator will illuminate if top plate temperature is over 50°C.
- For correct control media by the temperature probe, make sure the sensor to completely immerse into the media (more than 20mm deep). If the sensor is in the air due to evaporation of media, Errors will be shown on display.
- Probe operating temperature range is :
   B Class : Max 250°C.
   A Class : Max 400°C

**Configuration Mode** Various configuration settings can be adjusted under the configuration mode such as heating mode, timer mode, maximum heating rate, low and high temperature limits, temperature unit, temperature offsetting, temperature control mode, and so on.


You can enter this configuration mode by pressing the CONFIG Button when the unit is properly connected to a power outlet but not operating.

You can terminate the configuration mode by using one of the following methods:


Step 1: Press the RUN/STOP Button.



Step 2: Turn the stirrer knob to select YES and push it to confirm the termination as shown below.



Or turn the stirrer knob clockwise until you get the following ESC display and then push it to confirm the termination as shown below:



Note also that the configuration mode is automatically terminated if there is no button operation for more than 20 seconds.

Temperature Unit Conversion (°C ↔°F) You can select the desired temperature unit as follows: Step 1: Press the CONFIG Button



4-14


Section 4 Control Panel

# Temperature Unit Conversion (°C $\leftrightarrow$ °F)

Step 2: Select the temperature unit conversion by turning the Stirrer Knob as shown below.




Step 3: Select either Celsius (°C) or Fahrenheit (°F) by turning and pressing the Stirrer Knob.



Step 4: Save the selected temperature unit by using the Stirrer Knob as


shown below.

[When you need to save the changes]



If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO".

[When you don't need to save the changes]



Step 5: Terminate the configuration mode as described in Configuration Mode.


# Selection of the Timer Mode

Two different timer modes are provided: Immediate Timer Activation (ti1) and Delayed Timer Activation (ti2). In case of the Immediate Timer Activation, timer starts immediately after setting the timer. In case of the Delayed Timer Activation, on the other hand, timer is activated only when the set temperature is reached.


The timer mode selection procedures are as follows: Step 1: Press the CONFIG Button.




Step 2: Select the timer mode shown below by turning and pressing the Stirrer Knob.



Step 3: If you want to change the timer mode, then turn and press the Stirrer Knob.




Step 4: Save the desired timer mode by using the Stirrer Knob as shown below. [When you need to save the changes]



If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]

# Selection of the Timer Mode



Step 5: Terminate the configuration mode as described in Configuration Mode.

#### NOTICE

 Refer to Top plate's temperature influence on Heating bath and Erlenmeyer flask's temperature for understanding of temperature relation between top plate and samples.

Temperature Offsetting<br/>(OFST/OFSP\*)The temperature shown on the Actual Temperature Display is<br/>measured by a temperature sensor inside the unit. However, this<br/>temperature can be different from the temperature of the external<br/>temperature probe. If needed, you can offset such temperature<br/>differences at 0.1°C interval.

\* Offsetting of top plate is displayed as OFST and offsetting of the external temperature probe is displayed as OFSP.

#### NOTICE

• The allowed range of the temperature offsetting:

OFST (-10 ~ +50°C) / OFSP (-3 ~ +3°C)

- See the graph in Accessories which shows the temperature difference depending on the top plate temperature.
- Refer to Temperature Offsetting for the relation between top plate's actual temperature and the shown temperature on the display in case of the unload condition.

If, for example, the actual temperature of the top plate is 100°C but the displayed temperature is 95°C, you can match the displayed temperature with the actual temperature of the unit by selecting the offset value of +5°C and save it as described below:

Temperature Offsetting (OFST/OFSP\*)

Step 1: Press the CONFIG Button

un E

Step 2: Select the offsetting mode (OFST or OFSP) by turning the Stirrer Knob and press the Stirrer Knob as shown below.



[OFST- Top plate offsetting]

[OFSP- External temperature probe]

Step 3: Select the offset value by turning the Stirrer Knob.


The allowed ranges of the temperature offsetting are from -10 to  $+50^{\circ}$ C for OFST and from -3 to  $+3^{\circ}$ C for OFSP.



[OFST- Top plate offsetting]

[OFSP- External temperature probe]

Step 4: Save the desired offsetting value by pressing the Stirrer Knob as shown below. [When you need to save the changes]




[OFST- Top plate offsetting]


[OFSP- External temperature probe]

If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]

# Temperature Offsetting (OFST/OFSP\*)







[OFSP- External temperature probe]

Step 5: Terminate configuration mode as described in Configuration Mode.

# Maximum Heating Rate Setting (LIMt)


The higher the heating rate, the faster the heat up time but the wider the temperature overshoot and undershoot. Therefore, if you want to reduce the temperature fluctuation, you need to limit the heating rate to a certain degree by limiting maximum heating rate. [Same function as Temperature Setting].

If, for example, the current heating rate is 100% but you want to set the heating rate limit at 60%, do as follows:


Step 1: Press the CONFIG Button.



Step 2: Select Limit (LIMt) turning and pressing the Stirrer Knob.




Step 3: Change the limit to the desired value by turning the Stirrer Knob and press the Stirrer Knob.




## Maximum Heating Rate Setting (LIMt)

Step 4: Save the desired limit by using the Stirrer Knob as shown below. [When you need to save the changes]



If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]



Step 5: Terminate the configuration mode as described in Configuration Mode.


#### NOTICE

- When the heat up time is too slow, check whether the maximum heating rate is set too low.
- Maximum heating rate set point can only be applied in top plate's control mode, not the external temperature probe mode.

### Changing the High Temp Limits (t-H)

Note that the default setting of the high temperature limit [t-h] 350°C. If needed, however, you can set your own temperature limits. If, for example, you want to set the high temperature limit at 200°C, do as follows:

Step 1: Press the CONFIG Button.




Step 2: select the temperature limit mode shown below by turning the Stirrer Knob. Press the Stirrer Knob and check the current high temperature limit:


4-20

Section 4 Control Panel


# Changing the High Temp Limits (t-H)




Step 3: Change the limit to the desired value by turning the Stirrer Knob.



Step 4: Save the desired limit by using the Stirrer Knob as shown below. [When you need to save the changes]



If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]



Step 5: Terminate the configuration mode as describe in Configuration Mode.

#### NOTICE

- The allowed range of the temperature limit setting is 0°C ~ 350°C.
- After terminating the configuration modes, new temperature setting should be arranged within the range of the temperature limit.


# Changing the Low Temp Limits (t-L)

Note that the default setting of the low temperature limit [t-L] 0°C. If needed, however, you can set your own temperature limits. If, for example, you want to set the high temperature limit at 100°C, do as follows:


Step 1: Press the CONFIG Button.

ו רח בו
---------


Step 2: Select the temperature limit mode shown below by turning the Stirrer Knob. Press the Stirrer Knob and check the current low temperature limit:



Step 3: Change the limit to the desired value by turning the Stirrer Knob.




Step 4: Save the desired limit by using the Stirrer Knob as shown below. [When you need to save the changes]



If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]

Section 4 Control Panel

# Changing the Low Temp Limits (t-L)



Step 5: Terminate the configuration mode as described in Configuration Mode.

#### NOTICE


- The allowed range of the temperature limit setting is  $0^{\circ}$ C ~  $350^{\circ}$ C.
- However, low temperature limit [t-L] set point should not be more than high temperature limit [t-H] set point.

# Check Configuration Settings (LISt)

You can check the current configuration settings as follows: Step 1: Press the CONFIG Button.



Step 2: Select the configuration checking mode [LiSt] shown below by turning the Stirrer Knob.



Step 3: Press the Stirrer Knob and check the current settings by turning the Stirrer Knob.

The default configuration settings are as follows:


MODE	Default Setting	MODE	Default Setting
Unit	Cels	Time	Ti1
IFST/OFSP	0.0	Limt	100.0
t-H	350.0	t-L	0.0
Coef	Opti	Mmin	1

Step 4: After checking the settings, terminate the checking mode.

[Either by pressing RUN/STOP Button]:

[Or by pressing the Stirrer Knob]:

Check Configuration Settings (LISt) Step 5: Terminate configuration mode as described in Configuration Mode.




## Auto-tuning (PID Parameter Calibration)

The PID parameters for temperature control can be automatically tuned to take your specific operating circumstances into consideration.


If you select User Mode [Coefficient], the unit will be operated with the auto tuning value.


The auto-tuning (calibration) procedures are as follows:


Step 1: Press the CONFIG Button.



Step 2: Select the auto-tuning mode [CALi] shown below by turning the Stirrer Knob.







4-24

Section 4 Control Panel

# Auto-tuning (PID Parameter Calibration)


Step 4: Terminate configuration mode as described in Configuration Mode. To cancel the auto-tuning during the process, turn off the main power switch. If cancelled, the Optimal Mode will be selected automatically with an alarm.



Coefficient (COEF)

Five user-selectable temperature control modes are provided for your convenience.


Step 1: Press the CONFIG Button.



Step 2: Select temperature control mode selection [CoEF] shown at right by turning the Stirrer Knob.



Step 3: Turn the Stirrer Knob to check the five temperature control modes in order. And select the desired mode by using the Stirrer Knob as shown below.



## **Coefficient (COEF)**

[User-selectable temperature control modes]


	Providing the optimal balance between the heat up	
Optimal Mode	time to reach the target temperature and the	
	allowed fluctuation range of temperature overshoot	
	and undershoot (factory default)	
	Providing the fastest heat up time but the widest	
Fast Mode	fluctuation range of temperature overshoot and	
	undershoot	
	Providing the slowest heat up time but the	
Slow Mode narrowest fluctuation range of temperature		
	overshoot and undershoot	
	Allowing linear heating up to the target temperature	
	using the maximum heating rate and switching on	
Point Mode	and off the heater based on the target temperature	
	(Note that this mode shows the widest fluctuation	
	range of temperature overshoot and undershoot.)	
User Mode	Allowing the auto-tuned parameters to be used for	
	temperature control	

Step 4: Save the desired mode by using the Stirrer Knob as shown at right. [When you need to save the changes]



If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]

#### **Coefficient (COEF)**



Step 5: Terminate configuration mode as describe in Configuration Mode.

Default (DEFA)If you want to change the current settings back to the default settings<br/>described in Selection of Timer Mode, do as follows:


Step 1: Press the CONFIG Button.

โมกา
------


Step 2: Select the default setting [dEFA] shown below by turning the Stirrer Knob.



Step 3 : Press the Stirrer Knob and confirm the change to the default settings by using the Stirrer Knob as shown below: [When you need to save the changes]




If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes] Default (DEFA)




- Step 4: Terminate configuration mode as described in Configuration Mode.
- Motor Minimum (Mmin) You can limit the stirring speed to a certain point by setting the limit of maximum stirring speed. The range of speed adjustable is from 1 to 1000rpm.

Step 1: Press the CONFIG Button.




Step 2: Select the Motor Minimum [Mmin] shown below by turning the Stirrer Knob.




Step 3: Change the limit to the desired value by turning the Stirrer Knob

#### and press the Stirrer knob.




Step 4: Save it by using the Stirrer Knob as shown below.



#### Motor Minimum (Mmin)

If you don't want to save it, press the Stirrer Knob when display shows "SAVE NO". [When you don't need to save the changes]



Step 5: Terminate configuration mode as described in Configuration Mode.


Escape (ESC)

Terminate the configuration mode.

Step 1: Press the CONFIG Button.



Step 2: Select and press ESC by turning the Stirrer Knob.



# Section 5 Safety Device

(1) High and Low Temperature Limits

The high temperature limit of the top plate is set to 350°C to protect you and your media. However, you can set your own high and low temperature limits to reduce operation time and also to avoid inadvertent mistakes.

- (2) If the temperature of the top plate exceeds the set temperature range, the heater will be automatically turned on and off and the Heater LED will be turned on and off accordingly.
- (3) Multiple Overheat Prevention Measures (ERROR 1 / ERROR 3) Built-in overheat prevention circuit will turn off the heater if the plate Temperature reaches 450°Cfor any reason. In addition, if the temperature of the main body exceeds 85°C the overheat prevention circuit also stops heating to protect the motor and the PCB. You will be alerted by both audible and visible signals will be activated.
- (4) In such cases, turn off the power switch and disconnect the power cord first. Then allow your unit to cool down completely before operating it again.
- (5) Hot Top Warning Indicator

The top plate temperature can remain very hot for some time even after the heater is turned off. To prevent injury or fire under such circumstances, this equipment has a hot top warning indicator on the control panel. This indicator will illuminate if the top plate temperature is over 50°C. Even so, do not rely on this indicator alone for your safety.

(6) Temperature Probe Warning (ERROR6,7)

On temperature control mode if the external probe is in the air or is not enough immersed into the media, error message will be shown on display and heater will be automatically turned off.

(7) In such cases, turn off the power switch and disconnect the power cord first. Then allow the temperature probe to completely immerse into the media (more than 20mm deep) and operate it again.

#### NOTICE

• If any overheat prevention activates or warning sound is alarming, turn off the power switch and disconnects the power cord first. Then allow your unit to cool down completely before operating it again.

# Section 6 Maintenance

ltem	Inspection Interval	
item	Daily	Weekly
Connection status of power cord or plug	•	
Damages in power cord or plug	•	
Damages or cleanliness of top plate		
Cleanliness of main body and accessories		
Damages in switches, buttons, LED's, dial		
knobs	•	
Heating capability check (up to 350°C)	•	
Stirring capability check (up to 2000 rpm)	•	
Assembly status of all parts or accessories		•

#### **Cleaning Product**

#### 

- Never immerse this unit in water or any other liquid.
- Do not allow any liquid or wet material to get inside the unit when cleaning.
- Do not reconnect this unit to power outlets until all cleaned surfaces have dried.

#### 

- Do not use chlorine bleach, ammonia-based cleaners, abrasives, ammonia, or metal scouring pads when cleaning.
- During cleaning and general operation, take care not to scratch the surface of the ceramic-coated top plate as this could result in subsequent thermal breakage.

Cleaning Product	Always make sure to keep top plate, main body, and accessories clean. Dirt and other foreign substances can cause fire or electric shock. Before attempting cleaning,		
	<ul> <li>Disconnect the power cord from the power outlet and ensure that the instrument is cool enough,</li> <li>Wipe with a soft dry cloth first to remove any foreign matter and, if not enough,</li> <li>Wipe with a soft damp cloth or a sponge soaked in water or diluted neutral detergent when necessary.</li> </ul>		
	Note that cleaning is made much easier if spills are attended to promptly.		
Relocation	If you need to move the instrument to another place,		
	<ul><li>(1) Disconnect the power cord from the power outlet,</li><li>(2) Pack the instrument and its accessories into the original packaging or any other suitable container before moving.</li></ul>		
	▲ CAUTION • Pay attention to avoid mechanical shock or vibration while moving the instrument. Damages caused by mechanical shock or vibration may result in injury or fire.		
Keeping Product	If you know you will not use this unit for an extended period of time,		
	<ul><li>(1) Disconnect the power cord from the power outlet and</li><li>(2) Clean the instrument with soft cloth.</li><li>(3) Pack the instrument properly and make sure to store it in dry place.</li></ul>		

# Section 7 Trouble Shooting

Electrical Trouble	Causes	Solution	
	Unsuitable power supply	Meet the electrical requirements of this	
		instrument before use.	
		Find out why blackout or cut-off happened and	
	Power cut-off by a circuit	restore power. If there is a short circuit or	
	breaker or power blackout	leakage, trace the source of the problem and fix	
No power		it.	
		Reconnect the power cord firmly to the power	
	Loose power connection	outlet as well as to the power socket at the back	
		of the instrument.	
	Damages in power cords,	Replace the damaged part with a proper one	
	power outlets or plugs.		
	Internal circuit failure	Contact Fisher Scientific for service.	
		Disconnect all the appliances connected to the	
Repetitive tripping of the	Electrical overload	breaker first and reconnect them one by one to	
circuit breaker		find the reason for the overload.	
	Internal circuit failure	Contact Fisher Scientific for service.	
		If the temperature of the main body exceeds	
	Power cut-off by built-in	$85^\circ$ C, the built-in overheat prevention circuit	
No operation with power	overheat prevention	stops heating to protect the instrument.	
on	circuit	In such cases, let the instrument cool down for	
		some time before power reconnection.	
	Internal circuit failure	Contact Fisher Scientific for service.	

Trouble During Operation Causes		Solution
	Failure to push the	Push the RUN/STOP Button.
	RUN/STOP Button Set point value is lower	
	than present temperature value.	Adjust Set point lower than Present temperature value.
No heat	Power cut-off by built-in	Turn off the main power switch and wait until the
	overheat prevention	heater cools down. Then, turn on the main
	circuit	power switch.
	Button switch failure	Contact Fisher Scientific for service.
	Internal circuit failure	Contact Fisher Scientific for service.
	Probe is not enough	Soak probe into the media more than 20mm of
No. on to a plan	soaked into the media.	it.
No or too slow	Too low setting of the	Turn the heater knob clockwise to increase the
temperature change	heating level	heating level
during heating	Too much media	Reduce the media volume.
	Internal circuit failure	Contact Fisher Scientific for service.
	Too much media	Reduce the media volume or increase the rpm
		more gradually.
	High viscosity of the	Increase the rpm more gradually or replace the
Stir bar decoupling		stir bar with a new one with less friction
	media	resistance (e.g., cone type).
	Decreased magnetic strength of the stir bar	Replace the old stir bar with a new one.
Knocking noise during	Uneven bottom of the vessel	Use a vessel with thin and flat bottom.
stirring	Loosened internal parts	Contact Fisher Scientific for service.
Abnormal anond control	Too much media	Reduce the media volume.
Abnormal speed control operation	Internal circuit failure or damaged motor	Contact Fisher Scientific for service.
Knob malfunction	Damaged knob	Contact Fisher Scientific for service.
LED display malfunction Damage due to chemical spill or overheat		Contact Fisher Scientific for service.

#### Section 7 Trouble Shooting

Trouble During Operation Causes		Solution
	Probe connection failure.	Reconnect probe in right way.
Probe lamp is not ON	Damage due to chemical	Contact Fisher Scientific for convice
	spill or overheat.	Contact Fisher Scientific for service.

# Error Messages

If you see below "Error codes", turn off the power. Restart after having the unit rest enough.


Error Code	Cause		Effect
ERRO NO 1		If Top plate is	Heater off
ERRONOT		overheated	Healer OII
		If sensor detects a	
ERRO NO 2		sudden temperature	Heater off
	Related to	changes	
ERRO NO 3		If PCB is overheated	Heater off
ERRO NO 4	Heater system	If sensor is not	Heater off
EKKO NO 4		connected correctly	rieater on
		If temperature difference	
ERRO NO 5		has arisen between two	Heater off
		sensors	
ERRO NO 6		If temperature probe is	Heater off
	Related to	in the air	Tieater off
	temperature	Temperature probe is	
ERRO NO 7	probe	not soaked enough into	Heater off
		media	

NA Model	EU Model	Description
11676275	15336617	FS Temperature Probe (PT100, B Class, Max to 250°C)
11676278	15366617	FS Temp Probe Advncd (PT100, A Class, Max to 400 $^\circ$ C)
11676269	15366607	FS Heating Bath
11676270	15376607	FS Transparent Shield (PC)
11676271	15386607	FS Support Rod (12Ø support, 400mm, M10)
11676276	15346617	FS Clamp Holder (PP body, Ø12mm)
11676277	15356617	FS 3 Prong Clamp ( 60mm grip)

# Section 8 Accessories

# Assembly of Support Rod


Hand-screw the support rod(s) into the threaded hole(s) as shown in the diagram. To firmly tighten the rod, use a proper wrench.



Section 8 Accessories


# Transparent Shield & Heating Bath Assembly

When necessary, the transparent shield and/or the heating bath can be easily assembled as shown below:





**Transparent Shield** 





Heating Bath

	Item / MODEL	Isotemp RT AVCD HPS
4001/	US Catalogue Number	11676262
120V EU Catalogue Number		N/A
0001/	US Catalogue Number	11676261
230V	EU Catalogue Number	15326607
	Temperature Range (°C/°F)	Max. 350 / 662
	Heating plate, temperature range, Max. (°C/°F)	350 / 662
Llastar	Control	Feedback Control with PID
Heater	Function	Offset, Auto tuning
	Temperature display	up to 350°C (0.1°C resolution)
	Heater output, max. (W)	600
	Speed range (Guarantee rpm)	30~2,000
	Speed display resolution (rpm)	1
Stirrer	Stirring capacity, Max (L/ cu ft, H <sub>2</sub> O)	20 / 0.7
	Motor Type	BLDC
	Magnetic Stir Bar, Max. (Ø x L, mm/ inch)	8 x 40 / 0.31 x 1.57
Droho	Туре	PT100, B class
Probe	Accuracy (°C, at 0°C)	±0.3
Timer		2 type, 1min to 99 hr 59 min
Cofety	Hot top warning (°C/°F)	Warning lamp (50 / 122)
Safety	Overheat Prevention	Top plate, Main body, PCB
Dimensior	n of top plate (Ø, mm / inch)	140 / 5.5
Overall Dimension (W x D x H, mm/ inch)		161 x 290 x 100 / 6.34 x 11.42 x 3.94
Weight (kg/lbs)		2.8 / 6.17
Electric Requirements (230V,50/60 Hz)		3A
Electric Requirements (120V, 60 Hz)		5A
Maximum load (kg/lbs)		25 / 55.1
Matarial	Main Body	Aluminum
Material Top plate		Ceramic coated aluminum

# Section 9 Technical Specifications

X Unless otherwise specified, the above-mentioned data represent values at 25°C and 60% relative humidity.

X Fisher Scientific reserves the right to make changes in design and specification without prior notice.

#### **Disposing of the Unit**

Disposing of your instrument must be done in an environmentally responsible way if it has been potentially exposed to bio-agents or radioactive samples. Failure to follow stringent requirements for instrument disposal may lead to actions against you and your organization.

- (1) First, check with your laboratory or organization to ensure that you are following all the policies and procedures for disposal of laboratory equipments.
- (2) If not possible, contact your local governing body for regulations regarding disposal of laboratory equipments. It is highly recommended that you to find a local service provider that can properly dispose of your instrument.

9-2

## FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY

When used in laboratory conditions and according to these operation instructions and maintenance, this product is warranted for 24 months against defective materials or workmanship. The 24 month warranty period begins from the delivery date of this product.

For product quality or performance issues, contact Fisher Scientific Customer Service.

10-1

#### **North America**

United States 1-800-766-7000 fishersci.com Canada 1-800-234-7437 fishersci.ca

#### Europe

Austria: +43(0)800-20 88 40 at.fishersci.com

Belgium: +32(0)56 260 260 be.fishersci.com

Denmark: +45 70 27 99 20 fishersci.dk

Germany: +49(0)180 5258221 de.fishersci.com

Ireland: +959(0)1 885 5854 ie.fishersci.com

Italy: +39 02 950 59 478 it.fishersci.com

Finland +358(0)9 8027 6280 fishersci.fi

France: +33(0)388 67 14 14 fishersci.fr Netherlands: +31(0)20 4887 70 00 nl.fishersci.com

Norway: +47 22 95 59 59 fishersci.no

Portugal: +351 21 425 33 50 pt.fishersci.com

Spain: +34 002 239 303 es.fishersci.com

Sweden: +46 31-68 94 30 fishersci.se

Switzerland: +41(0)56 618 41 11 ch.fishersci.com

UK: +44(1)1509 555 500 fisher.co.uk

# Find out more at fishersci.com

© 2017 Thermo Fisher Scientific Inc. All rights reserved. Trademarks used are owned as indicated at fishersci.com/trademarks.

