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Series TM

User Instructions:

Standard Heating Mantles, Installation and Operation

With proper care and operation, Glas-Col mantles will give long and efficient service. Contamination, overheating and general misuse will greatly reduce the life of a mantle. The Glas-Col mantle consists of an insulated electrical-resistance heater and requires a controlled electrical source. The use of a variable voltage source such as our manual control or automatic controller is recommended to prevent overheating and to effect accurate temperature control. A selection of temperature and voltage controls are available from the Glas-Col web site.

The limitation in the operation of the Glas-Col mantle is the maximum temperature allowed for the heating zone. This is 400°C (752°F) for glass fabric in Series O, M, and TM series mantles; and 600°C (1112°F) for series STM mantles. Series O spherical and hemispherical mantles should be used with the appropriate Glas-Col support. (Cat. No. O-524 through O-550 extension supports or O-504 and O-505 tripod supports.) The use of these supports is the only method approved by Glas-Col, and we do not assume responsibility for mantles damaged as a result of being supported otherwise— for example, directly on ring stands, cork ring, lab jacks, etc.

For maximum safety, the three-prong plug of the cord on the Series M, TM and STM mantles should always be inserted into a grounded outlet. Since the exterior of the Series O mantles is an insulator, these mantle are furnished with a two-wire cords.

In 12-liter Series O spherical mantles, there are three heating circuits—two in the lower half and one in the upper half. The heat input in each circuit should be controlled with a suitable control device. The two lower circuits furnish heat for boiling the liquid contents, while the upper circuit prevents condensation of the vapors. It may often be necessary to operate the upper circuit at not more than 50% of full power. For low-boiling liquids, the upper circuit need not be used. When the flask is more than half full of liquid, the two bottom circuits may be operated at the rated voltage. When the liquid level falls below the halfway mark in the flask, the power to the upper circuit of the lower half should be reduced to 50% or less. This will prevent superheating of the vapors, and overheating of the glass fabric in this circuit. A similar procedure should be followed when using a Series M, TM or STM mantle with a Series MO heating top. The procedure is also applicable to other multiple-circuit mantles.

Operation: Glas-Col mantle operation requires **Manual** or **Automatic** temperature control to maintain the desired temperature and to prevent overheating of the mantle, the vessel and its contents. Low-temperature (type GF) mantles should not be operated above 400°C (752°F). High-temperature (type QF) mantles should not be operated above 600°C (1112°F). Automatic and manual controls are available from Glas-Col. Call the factory if assistance is needed in selecting the proper temperature control.

Maintenance: Glas-Col mantles do not require regularly scheduled maintenance. However, regular inspection is recommended, and damaged mantles should be removed from service immediately. Glas-Col mantles should be protected from chemical spillage, mechanical damage and corrosive atmospheres so far as possible. Damage from spillage in Series M or TM mantles can be reduced greatly by using the Glas-Col Splash Guard or Poncho Safety Shield.

Spare Parts: Spare parts for Glas-Col heating mantles are not generally available with the exception of fastener components, electrical connectors and cord sets. To stock other spare parts, a complete spare mantle is usually necessary.

Limited Warranty:

See the current Glas-Col warranty policy located under the General Sales Policy on the Glas-Col website at www.glascol.com

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Temperature Control
with Limit

Digital Indicating
Control

If your heating process requires basic control or several steps, Glas-Col has the control for you. We offer bench-top or custom designed industrial systems with features like Ramp/ Soak that allows up to a 40-step profile. You can profile Temperature, Time, Hold, Soak and End steps to create the ideal profile for your process. The built in adaptive control technology provides even tighter control for these demanding applications. Several input types are available. Visit our website for addition safety devices.

Glas-Col, LLC

711 Hulman Street, P.O. Box 2128, Terre Haute, IN 47802-0128

Phone: (812) 235-6167, Fax: (812) 234-6975

E-mail: pinnacle@glascol.com

Web: www.glascol.com