

EcoChyll X1 Benchtop Rotary Evaporator

Specifications:

Equipment

- Dimensions (approx.) – 12" L x 12" W x 29" H

Electrical

- Power cons. 645 watts
- Available in both 115 and 220 Volts, 50/60 Hz

Technical

- EcoChyll cooling capacity – *See brochure for full chart*
- EcoChyll minimum temperature – -40°C
- Metal parts – Stainless steel 316
- Glass parts – Borosilicate glass
- Condensation coils – coated with chemically-resistant fluoropolymer Halar
- Sealing and tubing materials – PTFE and FEP
- Warranty** – Two-Years Manufacturer (includes parts & labor, excludes all glass components)