

# Fisher Scientific Mini Vortex Mixer

*Mini in Size, Large in Power*



## Feature Overview

The Mini Vortex Mixer is an extremely compact vortex mixer with a horizontally circular, orbital mixing motion. The unique cup design accepts tubes and vessels of all sizes, including 15 and 50mL conical tubes.

The convenient touch function activates the motor immediately upon contact with a tube or vessel. Simply press the center cup head, and the mixer provides a smooth, instant vortexing of samples.

- Compact, powerful and reliable
- Convenient touch function
- Fixed Speed, 2800rpm
- Smooth, instant vortexing
- Quiet and stable



## Ordering Information

Description	Electrical	Cat. No.
<b>Mini Vortex Mixers</b>		
Mini Vortex Mixer	115V, 50-60Hz	14-955-151
Mini Vortex Mixer	230V, 50-60Hz	14-955-152

## Specifications


<b>Speed:</b>	Fixed, 2800rpm
<b>Orbit</b>	4.5mm
<b>Dimensions:</b>	3.8 x 3.8 x 2.5 in. 96 x 97 x 64mm
<b>Weight:</b>	1 lbs / 0.5kg
<b>Electrical:</b>	100 to 240V, 50-60Hz (0.6A)
<b>Warranty:</b>	2 Years

© 2013 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).

### In the United States:

For customer service, call 1-800-766-7000  
To fax an order, use 1-800-926-1166  
To order online: [www.fishersci.com](http://www.fishersci.com)

13\_0175 JJ/BC 10/13 BN0923134

 **Fisher  
Scientific**  
Part of Thermo Fisher Scientific