



Product description

Order number	62.547.100
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: red, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/rack

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Page 2

Size

Diameter	28 mm
Length of product	114 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	rack
Piece(s) / inner box	25
Piece (s) / outer case	200
Piece(s) / pallet	7200
Case volume	0.0436 cbm
Weight of product	0.0175 kg
Weight of case	4.01 kg
EAN of inner box	4038917376259
EAN case	4038917376242

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

